


Product Information

Disposable Temperature Probes

Probes for measurement of body temperature; for rectal, oesophageal, nasal and skin contact application




Contact us for detailed order information. OEM and private labelling on request.

Description

The disposable temperature probes are compatible to the YSI 400 series. They are high accuracy, easy-to-use probes with different plug types in sterile single package.

Technical data

Measuring principle:	NTC
Characteristic curve:	YSI 400 compatible
Measuring range:	0 - 50°C
Accuracy:	± 0,1°C in the range of 25 - 45°C
Standard Plug type:	Molex compatible female; other types on request



Sensatronic GmbH is a company certified according to ISO 13485 and Annex II of the Medical Device Directive 93/42/EEC (MDD).



Made in Germany 

Sensatronic GmbH
Am Ring 9
D-23970 Wismar, Germany
Tel. +49-3841-22 52 440-0
Fax +49-3841-22 52 442
info@sensatronic.com
www.sensatronic.com